



**PRODUCT RECALL**  
**PI-24-5089-FA-EXP**  
**Venclose Mixed Signal Wires**

18 July 2025,

Dear BD Distributor,

**Type of Field Action: Removal (Recall)**

**Affected Product**

This notice is limited to the catalog numbers / lot numbers listed in Attachment 1. No other catalog numbers or lot numbers are affected.

Product Name (Brand Name as per labelling)	Catalog No	Lot No., Serial No., or Software Version	UDI-DI	Expiration Date	Package Size
See Attachment 1					

Table 1: Affected Product Table

*Note: Lots highlighted Yellow have been added to the field action.*

**Description of the Problem:**

In October 2024, BD issued a Recall Notice for specific lots of Venclose™ RF Ablation Catheters, that may have been manufactured with internal wiring errors. This error in wire configuration would allow inaccurate temperature detection and may result in devices reaching temperatures higher than intended, while the user receives incorrect warnings that the catheter has not reached temperature.

Since the original recall notice was issued, BD received additional customer complaints, including complaints related to internal wiring issues and catheters arriving nonfunctional due to a lack of battery power that were not included in the original recall notice. BD is expanding the scope of this recall to include all lots of Venclose™ RF Ablation Catheters within expiration. Additionally, BD is updating the action from “Field Work” to a “Product Removal” for any remaining product inventory on the market.

**Clinical Risk Statement:**

The internal wiring error associated with the Venclose RF Ablation Catheter may result in excessive heat that could cause the following issues, potentially requiring additional medical care should they occur: heat-induced injury to nerves, skin burns, difficulty removing the device, damage to treated vessels or adjacent vessels, endothermal heat-induced thrombosis, deep vein thrombosis, pain, skin discoloration, or embolism of device components. Thrombosis and embolism of device components could lead to pulmonary embolism.



**Complaint & Adverse Event Statement:**

To date, eight complaints from health care providers related to the internal wiring error, three of which are related to a reportable serious injury event.

**Actions for Clinical Users:**

1. Cease use of the affected product and discard immediately.  
*Note: The catheter contains an embedded internal Alkaline battery that is not intended to be removed by the user. Dispose of the device, components and packaging materials using standard hospital practices and in compliance with all applicable regulations for electronics and batteries. Never incinerate discarded batteries.*
2. Possible adverse events associated with this issue are recognized complications from venous ablation procedures and patients are routinely monitored for these events. Patients who have undergone treatment with a Venclose™ RF Ablation Catheter that experienced difficulties in achieving the desired temperature may require additional follow-up beyond the usual standard of care at the discretion of the treating physician.

**Action To be Taken by BD:**

1. BD has placed product on hold and temporarily paused manufacturing of the Venclose™ RF Ablation Catheter. BD continues to investigate these issues and is committed to taking necessary actions to prevent recurrence.
2. BD will issue credit / product replacement upon receiving a completed customer response form.

**Product Distribution Time Frame:**

BD distributed the affected lots from 26 Sept 2022 through 4 June 2025. Our records indicate you may have received product from an affected lot.

**Please Take the Following Actions:**

1. Please discontinue distribution of product listed in this letter and immediately dispose of all affected product remaining in your possession in accordance with your local facility's process and in compliance with all applicable local, state, and federal regulations for electronics and batteries.
2. Identify all customers within your distribution network that purchased any affected product as defined in this notification. Provide a copy of the attached customer letter to all customers to advise them of this field action notification on BD's behalf.
3. Complete the attached Distributor Response Form and return to the BD contact noted on the form confirming acknowledgement of the recall notification, whether you have any of the impacted product, so that BD may acknowledge your receipt of this notification and process your credit.
4. Indicate on the response form the quantity of product identified at your facility and confirm that this inventory was destroyed.
5. Please contact your BD representative if you require assistance with this process or you may also contact the number below for further assistance.
6. Credit will be issued to distributors for all discarded products with complete Certificate of Destruction.
7. Report any adverse health consequences experienced with the use of this product to BD.



2 International Business Park Road  
The Strategy #08-08  
Singapore 609930  
Registration No. 201114149N  
[bd.com](http://bd.com)

BD is committed to advancing the world of health. Our primary objectives are patient and user safety and providing you with quality products. We apologize for any inconvenience this issue may have caused you and thank you in advance for helping us to resolve this matter as quickly and effectively as possible.

Yours Sincerely,

Signed by:

*Noorhafizoh Mdlazim*



Signer Name: Noorhafizoh Mdlazim  
Signing Reason: I approve this document  
Signing Time: 18-Jul-2025 | 12:21:15 PM SGT

662F2F310B9B4BB8BA5F81536E8D0F51

---

Noorhafizoh Mdlazim  
SEA Quality Lead



**DISTRIBUTOR RESPONSE FORM**  
**PI-24-5089-FA-EXP**  
**Venclose Mixed Signal Wires**

Please fill in the information below so that we may acknowledge your receipt of this notification. Complete and return the completed form to **SEA\_Quality** [SEA\\_Quality@bd.com](mailto:SEA_Quality@bd.com) / local BD representative by **25 July 2025**.

Please tick as appropriate.

- I have read and understood the attached notice taken appropriate actions
- We do not have affected product(s) in our inventory.
- We have identified the affected product(s) in our inventory. Affected products have been quarantined until disposal. Upon disposal, we will provide a copy of the Certificate of Destruction/proof of destruction to BD.

The expected date of destruction is: \_\_\_\_\_.

- We have identified all customers that purchased the affected catalog numbers and will notify the affected customers of this notice. The overview of the distribution to the customers are as attached in the Distribution Overview.

Product Description	Catalog No.	Lot No.	Quantity (eaches) received)	Quantity (eaches) sold	Remaining Quantity (eaches) in inventory to be *destroyed

\* Please provide a copy of the Destruction Certificate/ Proof of Destruction for the disposal of all affected units.



**Completed by:**

<b>Name:</b>	
<b>Signature:</b>	
<b>Date:</b>	
<b>Facility / Address / Telephone Number:</b>	



**Attachment 1: Affected Product Table**

Product Name	Catalog	Lot No.	UDI-DI	Expiration Date	Package Size
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66030181	(01)00858254006015(17)241006(10)66030181	10/06/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	65497423	(01)00858254006015(17)241009(10)65497423	10/09/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66599182	(01)00858254006015(17)241011(10)66599182	10/11/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66186400	(01)00858254006008(17)241013(10)66186400	10/13/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66634308	(01)00858254006015(17)241017(10)66634308	10/17/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66784768	(01)00858254006015(17)241020(10)66784768	10/20/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66802218	(01)00858254006008(17)241021(10)66802218	10/21/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66793696	(01)00858254006008(17)241024(10)66793696	10/24/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66793709	(01)00858254006015(17)241026(10)66793709	10/26/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66793710	(01)00858254006015(17)241028(10)66793710	10/28/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66863470	(01)00858254006008(17)241031(10)66863470	10/31/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66793697	(01)00858254006008(17)241102(10)66793697	11/02/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66793698	(01)00858254006008(17)241107(10)66793698	11/07/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	66981017	(01)00858254006008(17)241110(10)66981017	11/10/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67190839	(01)00858254006008(17)241115(10)67190839	11/15/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66897727	(01)00858254006015(17)241103(10)66897727	11/03/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	66947733	(01)00858254006015(17)241113(10)66947733	11/13/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	67054619	(01)00858254006015(17)241117(10)67054619	11/17/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67017611	(01)00858254006008(17)241120(10)67017611	11/20/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67241133	(01)00858254006008(17)241123(10)67241133	11/23/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67105253	(01)00858254006008(17)241125(10)67105253	11/25/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67429950	(01)00858254006008(17)241128(10)67429950	11/28/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	67308611	(01)00858254006015(17)241130(10)67308611	11/30/2024	Cases of 10 units



Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67429951	(01)00858254006008(17)241204(10)67429951	12/04/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	67429952	(01)00858254006008(17)241204(10)67429952	12/04/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	V67114001	(01)00858254006015(17)241108(10)V67114001	11/08/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	V67114002	(01)00858254006015(17)241108(10)V67114002	11/08/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	V67114003	(01)00858254006015(17)241108(10)V67114003	11/08/2024	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	69548248	(01)00858254006015(17)250322(10)69548248	03/22/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69548197	(01)00858254006008(17)250327(10)69548197	03/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69659895	(01)00858254006008(17)250331(10)69659895	03/31/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69548196	(01)00858254006008(17)250410(10)69548196	04/10/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69548198	(01)00858254006008(17)250416(10)69548198	04/16/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	69548249	(01)00858254006015(17)250419(10)69548249	04/19/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	69548250	(01)00858254006015(17)250423(10)69548250	04/23/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	69567357	(01)00858254006015(17)250427(10)69567357	04/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69659896	(01)00858254006008(17)250502(10)69659896	05/02/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69748757	(01)00858254006008(17)250505(10)69748757	05/05/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	69928277	(01)00858254006008(17)250509(10)69928277	05/09/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	69567358	(01)00858254006015(17)250511(10)69567358	05/11/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	70322501	(01)00858254006015(17)250515(10)70322501	05/15/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70487066	(01)00858254006008(17)250517(10)70487066	05/17/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70634659	(01)00858254006008(17)250519(10)70634659	05/19/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	70597994	(01)00858254006015(17)250523(10)70597994	05/23/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	70751594	(01)00858254006015(17)250525(10)70751594	05/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70673130	(01)00858254006008(17)250529(10)70673130	05/29/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70414634	(01)00858254006008(17)250601(10)70414634	06/01/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70868426	(01)00858254006008(17)250604(10)70868426	06/04/2025	Cases of 10 units



Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	70791660	(01)00858254006015(17)250606(10)70791660	06/06/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	70523716	(01)00858254006015(17)250607(10)70523716	06/07/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70933358	(01)00858254006008(17)250608(10)70933358	06/08/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	70994254	(01)00858254006008(17)250612(10)70994254	06/12/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71054985	(01)00858254006015(17)250613(10)71054985	06/13/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71075353	(01)00858254006008(17)250614(10)71075353	06/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71197363	(01)00858254006008(17)250618(10)71197363	06/18/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71197454	(01)00858254006015(17)250620(10)71197454	06/20/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71319092	(01)00858254006015(17)250622(10)71319092	06/22/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71278440	(01)00858254006008(17)250625(10)71278440	06/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71278441	(01)00858254006008(17)250628(10)71278441	06/28/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71361455	(01)00858254006015(17)250630(10)71361455	06/30/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71197455	(01)00858254006015(17)250703(10)71197455	07/03/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71197364	(01)00858254006008(17)250704(10)71197364	07/04/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71197365	(01)00858254006008(17)250705(10)71197365	07/05/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71609707	(01)00858254006008(17)250707(10)71609707	07/07/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71650329	(01)00858254006008(17)250710(10)71650329	07/10/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71589682	(01)00858254006015(17)250712(10)71589682	07/12/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71650328	(01)00858254006008(17)250714(10)71650328	07/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71818971	(01)00858254006008(17)250717(10)71818971	07/17/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71818972	(01)00858254006008(17)250719(10)71818972	07/19/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71780480	(01)00858254006015(17)250721(10)71780480	07/21/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71780481	(01)00858254006015(17)250724(10)71780481	07/24/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71629666	(01)00858254006008(17)250727(10)71629666	07/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71913951	(01)00858254006008(17)250731(10)71913951	07/31/2025	Cases of 10 units



Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71197574	(01)00858254006015(17)250801(10)71197574	08/01/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72167973	(01)00858254006015(17)250803(10)72167973	08/03/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	71197533	(01)00858254006008(17)250806(10)71197533	08/06/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72216758	(01)00858254006015(17)250809(10)72216758	08/09/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72282293	(01)00858254006015(17)250810(10)72282293	08/10/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72248924	(01)00858254006008(17)250807(10)72248924	08/07/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72151238	(01)00858254006008(17)250814(10)72151238	08/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72216513	(01)00858254006008(17)250816(10)72216513	08/16/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72216631	(01)00858254006015(17)250821(10)72216631	08/21/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72475111	(01)00858254006008(17)250823(10)72475111	08/23/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72405570	(01)00858254006008(17)250825(10)72405570	08/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72440542	(01)00858254006015(17)250828(10)72440542	08/28/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72576947	(01)00858254006008(17)250831(10)72576947	08/31/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72524414	(01)00858254006008(17)250903(10)72524414	09/03/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72609982	(01)00858254006015(17)250905(10)72609982	09/05/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72640519	(01)00858254006008(17)250907(10)72640519	09/07/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72768920	(01)00858254006008(17)250910(10)72768920	09/10/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72735944	(01)00858254006015(17)250912(10)72735944	09/12/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72491306	(01)00858254006015(17)250914(10)72491306	09/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72812308	(01)00858254006008(17)250918(10)72812308	09/18/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	72845915	(01)00858254006015(17)250920(10)72845915	09/20/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	71974718	(01)00858254006015(17)250821(10)71974718	08/21/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73048080	(01)00858254006008(17)250925(10)73048080	09/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73024457	(01)00858254006015(17)250927(10)73024457	09/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73184835	(01)00858254006015(17)250929(10)73184835	09/29/2025	Cases of 10 units



Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73141510	(01)00858254006008(17)250927(10)73141510	09/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73218082	(01)00858254006015(17)251005(10)73218082	10/05/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73310842	(01)00858254006244(17)251005(10)73310842	10/05/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73184825	(01)00858254006008(17)251008(10)73184825	10/08/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73319890	(01)00858254006015(17)251011(10)73319890	10/11/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73424182	(01)00858254006015(17)251012(10)73424182	10/12/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73441652	(01)00858254006008(17)251015(10)73441652	10/15/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73319925	(01)00858254006244(17)251009(10)73319925	10/09/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73493157	(01)00858254006008(17)251017(10)73493157	10/17/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73388088	(01)00858254006244(17)251023(10)73388088	10/23/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73542584	(01)00858254006008(17)251025(10)73542584	10/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	72978442	(01)00858254006008(17)250922(10)72978442	09/22/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	73319896	(01)00858254006251(17)251020(10)73319896	10/20/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73656671	(01)00858254006008(17)251029(10)73656671	10/29/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73441733	(01)00858254006015(17)251102(10)73441733	11/02/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	73573950	(01)00858254006251(17)251107(10)73573950	11/07/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73623196	(01)00858254006244(17)251109(10)73623196	11/09/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73923359	(01)00858254006244(17)251112(10)73923359	11/12/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73954408	(01)00858254006244(17)251114(10)73954408	11/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	73808592	(01)00858254006015(17)251121(10)73808592	11/21/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	73656672	(01)00858254006008(17)251123(10)73656672	11/23/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	74103378	(01)00858254006015(17)251125(10)74103378	11/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	74171115	(01)00858254006015(17)251127(10)74171115	11/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	74219123	(01)00858254006015(17)251129(10)74219123	11/29/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	74137006	(01)00858254006008(17)251201(10)74137006	12/01/2025	Cases of 10 units



Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	74137007	(01)00858254006008(17)251204(10)74137007	12/04/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	73890876	(01)00858254006251(17)251119(10)73890876	11/19/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73987507	(01)00858254006244(17)251116(10)73987507	11/16/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	74288518	(01)00858254006008(17)251206(10)74288518	12/06/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	74070677	(01)00858254006251(17)251214(10)74070677	12/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	74020834	(01)00858254006244(17)251214(10)74020834	12/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	73319926	(01)00858254006244(17)251214(10)73319926	12/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	74507761	(01)00858254006244(17)260107(10)74507761	01/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	74889457	(01)00858254006244(17)260110(10)74889457	01/10/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-100	74355356	(01)00858254006015(17)260112(10)74355356	01/12/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	74355347	(01)00858254006008(17)260114(10)74355347	01/14/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC-10A2.5-6F-60	74989575	(01)00858254006008(17)260121(10)74989575	01/21/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	74933204	(01)00858254006244(17)260116(10)74933204	01/16/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	75016007	(01)00858254006244(17)260124(10)75016007	01/24/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	75127795	(01)00858254006244(17)260130(10)75127795	01/30/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	75204432	(01)00858254006244(17)260204(10)75204432	02/04/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	75276852	(01)00858254006244(17)260206(10)75276852	02/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	74507759	(01)00858254006251(17)260208(10)74507759	02/08/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	75305699	(01)00858254006244(17)260214(10)75305699	02/14/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	75087737	(01)00858254006008(17)260215(10)75087737	02/15/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	74959653	(01)00858254006015(17)260222(10)74959653	02/22/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	75563807	(01)00858254006015(17)260227(10)75563807	02/27/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	75439758	(01)00858254006008(17)260304(10)75439758	03/04/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	75784764	(01)00858254006008(17)260306(10)75784764	03/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	75828822	(01)00858254006008(17)260311(10)75828822	03/11/2026	Cases of 10 units



Venclose® RF Ablation Catheter	VC10A256F60	75940087	(01)00858254006008(17)260314(10)75940087	03/14/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	76272657	(01)00858254006008(17)260326(10)76272657	03/26/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	76413778	(01)00858254006015(17)260402(10)76413778	04/02/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	76413837	(01)00858254006244(17)260404(10)76413837	04/04/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	76480490	(01)00858254006244(17)260406(10)76480490	04/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	76413830	(01)00858254006251(17)260409(10)76413830	04/09/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	76413711	(01)00858254006008(17)260412(10)76413711	04/12/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	76960626	(01)00858254006008(17)260422(10)76960626	04/22/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	77048331	(01)00858254006008(17)260424(10)77048331	04/24/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	77093851	(01)00858254006008(17)260425(10)77093851	04/25/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	77115614	(01)00858254006008(17)260502(10)77115614	05/02/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	76526332	(01)00858254006244(17)260506(10)76526332	05/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	77379776	(01)00858254006244(17)260507(10)77379776	05/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	77379777	(01)00858254006244(17)260507(10)77379777	05/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	76593887	(01)00858254006251(17)260507(10)76593887	05/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	77536423	(01)00858254006015(17)260524(10)77536423	05/24/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	77750390	(01)00858254006015(17)260529(10)77750390	05/29/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	77272719	(01)00858254006008(17)260603(10)77272719	06/03/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	77948128	(01)00858254006008(17)260605(10)77948128	06/05/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	77984934	(01)00858254006008(17)260611(10)77984934	06/11/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78094204	(01)00858254006008(17)260617(10)78094204	06/17/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78220237	(01)00858254006008(17)260621(10)78220237	06/21/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78295294	(01)00858254006008(17)260626(10)78295294	06/26/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78387756	(01)00858254006008(17)260701(10)78387756	07/01/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78481656	(01)00858254006008(17)260709(10)78481656	2026-07-09	Cases of 10 units



Venclose® RF Ablation Catheter	VC10A256F60	78634661	(01)00858254006008(17)260712(10)78634661	2026-07-12	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78692654	(01)00858254006008(17)260718(10)78692654	2026-07-18	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	77853714	(01)00858254006015(17)260722(10)77853714	2026-07-22	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78803219	(01)00858254006008(17)260729(10)78803219	2026-07-29	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78924909	(01)00858254006008(17)260801(10)78924909	2026-08-01	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78953690	(01)00858254006008(17)260806(10)78953690	2026-08-06	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79001814	(01)00858254006008(17)260813(10)79001814	2026-08-13	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	78857191	(01)00858254006015(17)260821(10)78857191	2026-08-21	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79073718	(01)00858254006008(17)260822(10)79073718	2026-08-22	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	79151080	(01)00858254006015(17)260827(10)79151080	2026-08-27	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79160974	(01)00858254006008(17)260828(10)79160974	2026-08-28	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	79208576	(01)00858254006015(17)260903(10)79208576	2026-09-03	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79217710	(01)00858254006008(17)260905(10)79217710	2026-09-05	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	79272760	(01)00858254006015(17)260910(10)79272760	2026-09-10	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	79345474	(01)00858254006015(17)260919(10)79345474	2026-09-19	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79364553	(01)00858254006008(17)260923(10)79364553	2026-09-23	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	79428613	(01)00858254006015(17)260925(10)79428613	2026-09-25	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	79479193	(01)00858254006015(17)260927(10)79479193	2026-09-27	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79462100	(01)00858254006008(17)261001(10)79462100	2026-10-01	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79531343	(01)00858254006008(17)261002(10)79531343	2026-10-02	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79539910	(01)00858254006008(17)261004(10)79539910	2026-10-04	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79557201	(01)00858254006008(17)261009(10)79557201	2026-10-09	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	78481657	(01)00858254006008(17)261014(10)78481657	2026-10-14	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79557202	(01)00858254006008(17)261017(10)79557202	2026-10-17	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79599378	(01)00858254006008(17)261025(10)79599378	2026-10-25	Cases of 10 units



Venclose® RF Ablation Catheter	VC10A256F60	79672390	(01)00858254006008(17)261029(10)79672390	2026-10-29	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79809444	(01)00858254006008(17)261104(10)79809444	2026-11-04	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79837661	(01)00858254006008(17)261107(10)79837661	2026-11-07	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79837662	(01)00858254006008(17)261126(10)79837662	2026-11-26	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79864789	(01)00858254006008(17)261126(10)79864789	2026-11-26	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79864790	(01)00858254006008(17)261126(10)79864790	2026-11-26	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	79837663	(01)00858254006008(17)261126(10)79837663	2026-11-26	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	80041650	(01)00858254006008(17)270102(10)80041650	2027-01-02	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	80041651	(01)00858254006008(17)270106(10)80041651	2027-01-06	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	80041652	(01)00858254006008(17)270115(10)80041652	2027-01-15	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	79359136	(01)00858254006251(17)260912(10)79359136	2026-09-12	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	79373868	(01)00858254006244(17)260616(10)79373868	2026-06-16	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	79762295	(01)00858254006244(17)261115(10)79762295	2026-11-15	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	79934996	(01)00858254006244(17)261121(10)79934996	2026-11-21	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	79987085	(01)00858254006244(17)261126(10)79987085	2026-11-26	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	79934993	(01)00858254006251(17)261218(10)79934993	2026-12-18	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	80689599	(01)00858254006244(17)270123(10)80689599	01/23/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	80689579	(01)00858254006008(17)270224(10)80689579	02/24/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	80572371	(01)00858254006008(17)270228(10)80572371	02/28/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81317621	(01)00858254006008(17)270303(10)81317621	03/03/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81317622	(01)00858254006008(17)270307(10)81317622	03/07/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81441841	(01)00858254006008(17)270311(10)81441841	03/11/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81380595	(01)00858254006008(17)270317(10)81380595	03/17/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	80794909	(01)00858254006244(17)270319(10)80794909	03/19/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81536712	(01)00858254006008(17)270321(10)81536712	03/21/2027	Cases of 10 units



Venclose® RF Ablation Catheter	VC10A256F100	80256667	(01)00858254006015(17)270327(10)80256667	03/27/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	81879617	(01)00858254006244(17)270331(10)81879617	03/31/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81923911	(01)00858254006008(17)270404(10)81923911	04/04/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	82246651	(01)00858254006008(17)270409(10)82246651	04/09/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	82356624	(01)00858254006008(17)270415(10)82356624	04/15/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	82487962	(01)00858254006008(17)270421(10)82487962	04/21/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F60	80321036	(01)00858254006251(17)270424(10)80321036	04/24/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	82621861	(01)00858254006008(17)270428(10)82621861	04/28/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	81923912	(01)00858254006008(17)270428(10)81923912	04/28/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	82155261	(01)00858254006244(17)270506(10)82155261	05/06/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	80321033	(01)00858254006015(17)270513(10)80321033	05/13/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	82062350	(01)00858254006015(17)270516(10)82062350	05/16/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VCOUS6F100	82965433	(01)00858254006244(17)270520(10)82965433	05/20/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	82768336	(01)00858254006008(17)270523(10)82768336	05/23/2027	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	073808592	(01)00858254006015(17)251121(10)073808592	11/21/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	073024457	(01)00858254006015(17)250927(10)073024457	09/27/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	073218082	(01)00858254006015(17)251005(10)073218082	10/05/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	074103378	(01)00858254006015(17)251125(10)074103378	11/25/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	075563807	(01)00858254006015(17)260227(10)075563807	02/27/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	074355356	(01)00858254006015(17)260112(10)074355356	01/12/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	076413778	(01)00858254006015(17)260402(10)076413778	04/02/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	074070677	(01)00858254006251(17)251214(10)074070677	12/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	076413830	(01)00858254006251(17)260409(10)076413830	04/09/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	075276852	(01)00858254006244(17)260206(10)075276852	02/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	076593887	(01)00858254006251(17)260507(10)076593887	05/07/2026	Cases of 10 units



Venclose® RF Ablation Catheter	VC10A256F100	074020834	(01)00858254006244(17)251214(10)074020834	12/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	074507759	(01)00858254006251(17)260208(10)074507759	02/08/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	075204432	(01)00858254006244(17)260204(10)075204432	02/04/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	076413837	(01)00858254006015(17)260404(10)076413837	04/04/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	076480490	(01)00858254006244(17)260406(10)076480490	04/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	077379776	(01)00858254006244(17)260507(10)077379776	05/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	076526332	(01)00858254006244(17)260506(10)076526332	05/06/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	077379777	(01)00858254006244(17)260507(10)077379777	05/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	074889457	(01)00858254006244(17)260110(10)074889457	01/10/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	075127795	(01)00858254006244(17)260130(10)075127795	01/30/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	074933204	(01)00858254006244(17)260116(10)074933204	01/16/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	075016007	(01)00858254006244(17)260124(10)075016007	01/24/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	075305699	(01)00858254006244(17)260214(10)075305699	02/14/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	073319926	(01)00858254006244(17)251214(10)073319926	12/14/2025	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	074507761	(01)00858254006244(17)260107(10)074507761	01/07/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F100	074959653	(01)00858254006015(17)260222(10)074959653	02/22/2026	Cases of 10 units
Venclose® RF Ablation Catheter	VC10A256F60	073890876	(01)00858254006251(17)251119(10)073890876	11/19/2025	Cases of 10 units